

EMR-validated column scale with eye-level display and Wi-Fi function

seca
Precision for health

seca 797

- + Medical-grade, robust scale for years of use with high patient volume
- + Digital measuring rod with an extra long measuring range
- + Two-line display shows height and weight
- + Oversized, low-profile full metal platform
- + 8-year warranty (4x the industry standard)
- + Optional: handrail seca 477 provides extra patient security
- + Optional: panda bear facilitates measuring and weighing of children
- + EMR-validated with built-in Wi-Fi module



EMR

Image contains barcode scanner from a third-party vendor.

EMR-validated column scale with eye-level display and Wi-Fi function



The sophisticated way of integrating weight and height

With its built-in Wi-Fi module the seca 797 is the perfect fit for seamless network integration of weight and height into your EMR (electronic medical record). In addition to information system connectivity, the seca 797 boasts a tilt-proof full metal platform and a stable single piece column with a display at eye-level. It is a quality every patient will notice, as soon they set foot on the platform. The scale is easy to move by staff due to the large rubber coated castors and a wide wheelbase.

Scale

Capacity	550 lbs / 250 kg
Graduation	0.2 lb / 100 g
Measuring range	4 – 90 in / 10 – 230 cm
Graduation	1/8 in / 1 mm
Dimensions (WxHxD)	16.1 x 53.4 x 22.7 in / 410 x 1,356 x 576 mm
Platform (WxHxD)	13.8 x 2.3 x 16.5 in / 350 x 59 x 418 mm
Net weight	28.9 lbs / 13.1 kg
Power supply	Batteries (optional), Power adapter
Functions	TARE, HOLD, Auto-HOLD, lbs/kg switch-over, lbs/kg switch-over lock, cm/inch switch-over, Adjustable damping, Auto-BMI, Wi-Fi
Interfaces	USB for barcode scanner, Wi-Fi
Accessories	+ Panda seca 459 + Handrail seca 477

Barcode scanner not included and may be purchased from a third-party vendor.

Measuring rod (included)

Measuring range (total)	4 – 90 in / 10 – 230 cm
Measuring range (digital)	48 – 90 in / 122 – 230 cm
Graduation	1/8 in / 1 mm



Optional handrail for increased patient stability

To accommodate patients in need of extra support during the weighing process, the seca 797 is equipped with the option to add an integrated handrail (seca 477). The easy-to-mount handrail raises patient security to a new level. The digital measuring rod is firmly attached to the column and provides a wide measuring range that is suitable for pediatrics to tall individuals. The two-lined display shows weight and height simultaneously and automatically calculated patients' BMI. The damping system and Auto-HOLD function enable a quick and precise measurement.



seca integration solutions: How workflows can clearly be improved.

The seca 797 can be integrated into your information system via the network. At the push of a button, measured values can be sent directly to the EMR, where they are immediately available to all users. Extended with USB barcode scanner (available from third-party vendor), both patient and user are identified and the measurement results are always correctly assigned. This saves time, increases process reliability and eliminates transmission errors.

Further information on seca integration solutions at:
integration.seca.com

Handrail seca 477 (optional)

Dimensions (WxHxD)	16.9 x 44.2 x 12.2 in / 429 x 1,122 x 309 mm
Net weight	8.8 lbs / 4 kg